

## 1. Preliminary remarks



Please read these instructions for use carefully before first use as correct temperature measurement depends only on the appropriate use of the device. These instructions describe the individual steps of temperature measurement using Infrared Thermometer and contain important and helpful hints for the reliable determination of the body temperature. Be sure to keep these instructions for use for future reference.

## 2. Introduction

Dear Customer, we are pleased that you have decided to purchase a clinical thermometer. Non-contact Infrared Thermometer is a quality product for measuring the human body temperature on the forehead. It is ideally suited for measurements on children as from 6 months of age, but the thermometer can be used also on adults. Given correct application the device ensures a fast and precise measurement of the body temperature in a very comfortable manner. We wish you all the best for your health. Indications for Use: Infrared Thermometer is intended to detect body temperature from forehead in the population including infant (above 6 months), child, adolescent, and adult.

## 3. BENEFITS OF THE THERMOMETER

### Non-contact measurement in a few seconds

High-tech innovative infrared sensor: to complete measurement safely and healthily in a few seconds.

### Body and Object temperature measurement

You can choose to perform measurement of either body temperature or object temperature by moving the slide switch.

### Fever Prompt

Red backlight display and 10 short beeps to warn the patient that he/she may have a fever.

### Display 30 measurement readings

When you enter the memory mode, you can read the last 30 measurement readings.

### Automatically display memories

This thermometer will display the last reading automatically for two seconds when starting the device.

## 3. IMPORTANT SAFETY INSTRUCTIONS

● Never use the product for other purposes except its original use please comply with the general safety precautions when used for children.

● This thermometer should never be immersed in water or other liquid (not waterproof). For cleaning and storage please read carefully the instructions in section 8.

● Please store the thermometer in a clean, dry environment; avoid direct sunlight; storage temperature should be between -4 to 131 / -20 to 55°C.

● Please do not use the thermometer if its sensing head or body have been damaged, and do not attempt to repair it when damaged.

● This non-contact thermometer is made up of high quality precision part. Do not drop the instrument! Protect it from severe impact and shock. Please do not twist its body and sensing head.


● Properly dispose of batteries, keeping them from small children and heat.

● If batteries are swallowed consult a doctor immediately.

### Warning:

● Using the Non-contact Thermometer should not replace visiting a doctor.

## 5. DESCRIPTION OF PRODUCT



No	Item	Execution status of mode	Description
1	Display screen		Display measurement readings And other corresponding symbols
2	SCAN button	Off status	● Press once to enter memory mode
		On status	● Press to take measurement
3	O/I button		Press to turn on/off thermometer
4	Battery cover		Protect the battery
5	Sensing head		Infrared sensing for measurement
6	Slide switch	Switch Measurement mode	Selected body temperature Measurement mode by sliding the switch upwards. Selected object temperature measurement mode by sliding the switch downwards.

### Compenets:



2"AAA"Batteries



Instruction Manual

## 6. DISPLAY SCREEN AND IMAGE DESCRIPTION

### Description of all kinds of images in display area

Image	Setting	Description
	Body mode	Can perform measurement of body temperature
	Object mode	Can perform measurement of object temperature
	Fahrenheit degree	Measure by degree centigrade
	Degree centigrade	
	Memory mode	Display measurement value from memory

### Different types of information in display area

Display screen	Procedure	Description
	Starting up	Turn on this device by Pressing O/I button. All information will be shown for two seconds.
	Memory	The last reading will automatically display on the screen for two seconds, containing M image, measurement mode image (human body) or object.
	Preparation of measurement	This device has been prepared for measuring, image of °F (°C) will keep twinkling.
	Finish of measurement	Reading will be displayed on the LCD screen, red backlighting for measurement value ≥ fever point, green backlighting for measurement value < fever point.
	To enter the next measurement	Around two seconds later, the image of °F (°C) will twinkle, this device is ready for the next measurement.


## 7. USING THE THERMOMETER

### How to Switch °C and °F?

Press the SCAN button about 8sec when the machine in off status, then it will occur °C or °F, if this unit not what you want, please press the SCAN button again immediately.

### Measurement of human body temperature

- Slide the switch upwards to enter body temperature measuring mode.
- Press O/I button, all images on screen will display.
- The last reading and M image will automatically display for two seconds.
- The thermometer will enter the ready state:
  - A short beep.
  - Image of °F (°C) and reference position will twinkle on the screen.



- To start the measurement, aim the thermometer at the center of the forehead no more than 5 cm away. If there is hair, sweat or dust on the forehead, please remove and wipe clean in order to improve accuracy.


- Press the "SCAN" button for two seconds, release and the measurement is started.

- You will hear a long beep, that means the measurement has finished, you can read the temperature value on the screen, as shown below.

- Intelligent analysis of body temperature.
  - If the measured temperature is below 99.5°F (37.5°C), the measurement result will display on screen with a long beep and green backlight.
  - If the measured temperature is higher than 99.5°F (37.5°C), the measurement result will display on screen with 10 short beeps. The back light display and the short beeps alert the patient that he/she has a temperature above 37.5°C.

### Measurement of the object temperature

- Slide down the switch to enter object temperature measuring mode.
- Press O/I button, all images on screen will display.
- The last reading and M image will automatically display for two seconds.
- The thermometer entered the ready state:
  - A short beep.
  - Image of °F (°C) and reference position will twinkle on the screen.



- Aim the thermometer at the center of object no more than 5 cm away. If there is water, dust or dirt on the object, please wipe clean in order to improve accuracy.

- Press the "SCAN" button for 2 seconds, while aiming at the center of the object, release and the measurement started.
- You will hear a long beep, that means the measurement has finished, you can read the temperature value on the screen, as shown below.

### Attention:

- Please keep users and the thermometer in a stable room temperature for at least 30 minutes before use.
- After continuous measurement, please wait at least two minutes to turn on or off the thermometer again.
- Do not take a measurement immediately or a period of time after nursing.
- Never use the thermometer in a high temperature environment.
- Before or during measurement, do not drink, eat, or move.
- Please clean the scanning area and remove dust, hair or sweat before using the thermometer.
- Please remove water, dust or dirt on the object before measurement.

- Please do not remove the thermometer before hearing the long beep at the end of measurement.

- If the sensing head is found to be dirty after or during measurement or before putting the thermometer back into the box, use alcohol cotton to wipe it carefully please.

- Please try to take measurement in the same area; results will be different in different areas.

- Because the temperature regulating function of baby is not perfect, do not take measurement immediately after the baby go into the environment with large temperature difference, to avoid temperature measurement deviation.

- Suggest using general type thermometers for measuring in the following cases:

- Measurement results are higher or lower than expected.
- Newborn baby younger than 100 days.
- Children under 3 years old who have a defective immune system, have very serious fever or no fever phenomenon.

## 8. USING THE THERMOMETER

The thermometer can automatically memories 30 measurement values. Please follow the following sequence to read memory value.

- In the turn-off state, press the "SCAN" button about 3 sec, you can enter the reading memory mode. The image "M" will flash on the screen.

- Press the "SCAN" button again, "1" will appear and then the last measured value and "M" image will flash.

- After that, press the "SCAN" button again, "2" will appear and then the penultimate measurement value and "M" image will flash.

- Continuously pressing "SCAN" button, 30 measurement values will be displayed in order.

- 30

- 30

Press the measurement button "SCAN" a gain after reading the thirtieth measurement value, it will return to the most recent measurement value.

## 9.ERROR MESSAGE

The screen display	Significance of display	Possible causes and Debugging method
	The ambient temperature is too high	The measuring temperature is higher than: ● 109.2°F (42.9°C)(body mode) ● 212°F (100°C)(object mode)
	The ambient temperature is too low	The measuring temperature is below: ● 89.6°F (32°C) ( body mode) ● 32°F (0°C) ( object mode)
	Measurement temperature is too high	The environment temperature is higher than the 104°F (40.0°C)
	Measurement temperature is too low	The environment temperature is below: ● 59°F (15°C)(Body mode) ● 41°F (5°C)(Object mode)
	Error display	System failure
	Blank screen	Please check whether the battery is properly installed and to check the negative and positive and positive poles of the battery
	No battery indicator	If the display screen shows only a fixed battery image, the battery should be replaced

## 10.CLEANING AND STORAGE

Clean the shell and sensing head with a cotton cloth wetted or moisten with the 70% isopropyl alcohol and make sure not to let the liquid into the inner part of the product. We advise that you should clean your thermometer when you finished your personal measurements every time.

Do not use abrasive cleaning agent, diluent or gas to clean. Do not immerse the product in water or other liquid.

Be careful not to scratch the LCD surface. If not used for a long time, please remove the battery from the apparatus in order to avoid damaging the thermometer from battery leakage.




## 11.Power supply and disposal notes

- Use only high-quality batteries (see specification in Chapter "Technical data").
- Never mix old and new batteries or batteries made by different manufacturers.
- Immediately remove exhausted batteries.
- If you do not intend to use the device for a longer period of time, you should remove to prevent possible leakage.

● Please help protect the environment: Do not dispose of batteries in household waste! Designated collection points or municipal collection / recycling centers to dispose of waste batteries.

This symbol on products and / or accompanying documents means that products must not be mixed with regular household waste.

This product comes with two 1.5V AAA batteries. When the LCD displays the battery symbol "▼", it is time to replace the batteries. As shown in the following picture, slide the battery cover off the end of the thermometer, insert new equivalent battery type and carefully replace the battery cover.



## 12. TECHNICAL SPECIFICATIONS

Product name: Non-contact Thermometer

Accuracy: Body mode:  $\pm 0.36^{\circ}\text{F}$  ( $95^{\circ}\text{F} \sim 107.6^{\circ}\text{F}$ ) /  $\pm 0.2^{\circ}\text{C}$  (35.0~42.0°C),  $\pm 0.54^{\circ}\text{F}$  ( $95^{\circ}\text{F} \sim 107.6^{\circ}\text{F}$  Outside) /  $\pm 0.3^{\circ}\text{C}$  (35.0~42.0°C Outside)

Object mode:  $\pm 0.9^{\circ}\text{F}$  ( $32^{\circ}\text{F} \sim 212^{\circ}\text{F}$ ) /  $\pm 0.5^{\circ}\text{C}$  (0~100°C)

Measuring range: Body mode: 89.6°F to 109.22°F / 32.0°C to 42.9°C

Object mode: 32°F ~ 212°F / 0°C to 100°C

Operating temperature: Body mode: 59°F to 104°F / 15°C to 40°C

Relative humidity: 20%~85%RH

Object mode: 41°F to 104°F / 5°C to 40°C

Screen: Liquid crystal display screen, display unit 0.1" / 0.1" C.

- Sound:
- Device start and preliminary measurements: 1 short "beep" sound.
  - Complete the measurements: one long "beep" sound.
  - System error or fault: three short "beep"
  - Measurement process: fast and slow "beep" sound.

- Memory:
- Automatically displays the last measurement reading.
  - Memory stores 30 readings.

Backlight display:

- When the device starts up, green backlight flashes on the screen for four seconds,
- When the measurement is completed, the green backlight will flash for five seconds for measurement values below 99.5°F / 37.5°C.
- When the measurement is completed, the red backlight will flash for five seconds for measurement values above 99.5°F / 37.5°C.

Storage/transportation temperature: -4°F to 131°F / -20°C to 55°C

Relative humidity: <85%RH

Automatic shutdown: After about three minutes without operation.

Battery: DC 3V (DC 1.5V AAA 2PCS)

Size: 153mm (L) X 36mm (W) X 36.8mm (H)

Weight: About 80g (with the battery)

## 13. WARRANTY

We grant you 2 years warranty after date of purchase. This product has been produced with the greatest care according to international quality standards, established in the European Guideline for Medical Products 93/42/EEC. The unit satisfies the requirements of ISO 80601-2-56 and ASTM E1965 and it was subjected to strict testing before delivery. Should you nevertheless have reasons of complaint, please send the thermometer together with the warranty card, filled out, to the service address given on the back.

Any damage caused by improper handling shall not be covered by the warranty. Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted within the warranty period. Be sure to include: date of purchase, dealer stamp, and name and address of responsible dealer.

## 14. EMC Declaration

Guidance and manufacturer's declaration-electromagnetic immunity.

The FR series device is intended for use in the electromagnetic environment specified below. The customer or the user of the FR series device should assure that they are used in such environment.

Immunity Tests	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge(ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	$\pm 8$ kV contact $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U <sub>T</sub> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U <sub>T</sub> ; 1 cycle and 70% U <sub>T</sub> ; 25/30 cycles Single phase: at 0° 0% U <sub>T</sub> ; 250/300 cycle	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VST300 device requires continued operation during power mains interruptions, it is recommended that the VST300 device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30A/m Magnetic field IEC 61000-4-8	30A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

NOTE: U<sub>T</sub> is the a.c. mains voltage prior to application of the test level.

## 15. EMC Declaration(Continued)

Guidance and manufacturer's declaration-electromagnetic immunity.

The FR850 series device is intended for use in the electromagnetic environment specified below. The customer or the user of the FR850 series device should assure that they are used in such environment.

Immunity Tests	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted FR IEC 61000-4-6	3 V rms 150 kHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz	N/A	Portable and mobile FR communications equipment should be used no closer to any part of the FR series device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter <b>Recommended separation distance</b> $d = 12 \sqrt{P}$
Radiated FR IEC 61000-4-3	10V/m 80MHz to 2.7GHz	10 V/m 80 MHz to 2.7 GHz	$d = 12 \sqrt{P}$ 80MHz to 800MHz $d = 23 \sqrt{P}$ 800MHz to 2.5GHz Where p is the maximum output power rating of the transmitter in watts (W)

			according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed FR transmitters, as determined by an electromagnetic site survey: " should be less than the compliance level in each frequency range." Interference may occur in the vicinity of equipment marked with the following symbol: 
--	--	--	---

NOTE1 At 80MHz and 800MHz, the higher frequency range applies.  
NOTE2 These guidelines may not apply in all situations.  
Electromagnetic is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VST300 device are used exceeds the applicable RF compliance level above, the VST300 device should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the VST300 device.

b) Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m

## 16. EMC Declaration(Continued)

Recommend separation distance between portable and mobile FR850 communications equipment and the FR850 series device

The FR850 series device is intended for use in an electromagnetic environment in which radiated FR850 disturbances are controlled. The customer or the user of the FR850 series device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile FR850 communications equipment (transmitters) and the FR850 series device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter w	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80MHz to 800MHz $d = 1.2 \sqrt{P}$	800MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	/	0.12	0.23
0.1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters(m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer.

NOTE1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# JIACOM

## Non-contact Thermometer

● FR850

## Instructions for use of products



Users must read the instructions carefully before using this product

## 17. Disposal

Once the product life has ended or his components do not work anymore their disposal should be carried out according to the current regulations.

## 18. Symbols or Abbreviation

- Type BF applied parts
- Manufacturer
- ATTENTION
- Follow instructions for use
- Serial number
- European Authorized Representative

**CE 0197** Notified Body

**EC REP** Shanghai International Holding Corp.GmbH(Europe)  
Eiffelstrasse 80, 20537 Hamburg, Germany

Shenzhen JIACOM Technology Co., Ltd  
Add: 301, No. 596-4 Dahe Village, Guancheng Community, Guanhu Street, Longhua District, Shenzhen, Guangdong, CHINA



DocID:FRIFU-01  
Version number:V01